

Part no.	311-8255200-060
Product name	機器說明書/MbH/英/.標準品/
Spec	L288*W250mm/雙面/雜誌紙/6折 (短邊彈簧3折+長邊對3折)完成尺寸L36*W62.5mm/65P/黑
Designer	Emily
Color	K100 K20

TROUBLESHOOTING		
Symptom	Possible Causes	Solutions
Fail to turn on the oximeter.	The batteries are dead.	Replace all batteries.
	The batteries are installed incorrectly.	Verify battery orientations.
SpO ₂ or pulse rate displays are missing.	Defective LCD displays.	Displayed values may not be reliable; discontinue use of the oximeter.
SpO ₂ or pulse rate displays unstably.	The finger might not be steady or placed incorrectly on the probe.	The finger might not be steady or placed incorrectly on the probe.
Disruption in the oximeter performance.	Electromagnetic interference (EMI).	Remove the oximeter from the EMI environment.
Battery is low and "  bAt Lo " is on display.	The batteries are low.	Replace the batteries immediately.
Backlight turns to blinking red (visual alarm is activated).	Oxygen saturation value is below 85%.	Consult healthcare professional immediately.

SPECIFICATION

Model No.: TD-8255
Dimension & Weight: 63(H) x 37(W) x 32(D) mm, 40 g without batteries
Display: LCD
Power Source: Two 1.5V AAA alkaline batteries
Battery Life: Batteries can be used continuously for 8 hours (for reference only. It depends on different brands of AAA alkaline batteries.)
External Output: Bluetooth
Measurement and Displayed Range: 0% to 100%
Increment: 1%
Accuracy: 100% to 80%: ±2%; 79% to 70%: ±3%; others are undefined.
Method: Dual wavelength LED
Operating Conditions: 5°C to 40°C (41°F to 104°F); 15% to 93% R.H. (non-condensing)
Storage / Transportation Conditions: -25°C to 70°C (-13°F to 158°F); 15% to 93% R.H. (non-condensing)
Product Warranty: At least 1 year (warranty period depends on local jurisdiction)
Range of Peak Wavelengths: 660 nm and 880 nm
Maximum Optical Output Power of Light Emitted by Oximeter Probe: 100 mW
Mode of Operation: Spot Check / Monitoring
Pulse Rate
Measurement and Displayed Range: 30 to 250 beats per minute (bpm)
Resolution: 1 bpm
Accuracy: ±1 bpm or ±1%, whichever is greater
Classification
Degree of Protection: Type BF Applied part
Safety: IEC60601-1
EMC: IEC60601-1-2
Harmonized Standard: ISO 80601-2-61
Home Use Medical Device: IEC60601-1-11
Water-resistance: IP22
Mode of Operation: Spot check / Monitoring

SYMBOL INFORMATION

	Serial number
	Alarm, general
	CE mark
	Unique device identifier
	This device does not belong to household waste and must be returned to a collection point for recycling electric and electronic devices according to local laws. If it contains batteries, the batteries should be removed and disposed in accordance with local regulations for separate collection of spent batteries.


CLINICAL PERFORMANCE

Tables below show Arms values measured using MbH TD-8255 Fingertip Pulse Oximeter in a clinical study. The individual and pooled measured Arms values in the discrete SpO₂ ranges of all 14 subjects are reported.

Subject	70% - 80% SaO ₂		80% - 90% SaO ₂		90% - 100% SaO ₂	
	Mean Bias	Arms	Mean Bias	Arms	Mean Bias	Arms
1	-1.00	1.89	1.25	1.80	0.00	1.03
2	1.27	1.71	-0.17	1.58	-0.81	1.20
3	-2.00	2.00	-1.90	1.97	-1.15	1.33
4	2.17	2.27	1.14	1.51	0.81	1.64
5	-1.11	2.11	1.25	1.94	-0.74	1.26
6	-0.57	2.07	-1.25	1.94	-1.00	1.22
7	1.00	1.78	2.00	2.00	0.20	1.06
8	1.30	1.97	0.50	0.71	-0.64	1.16
9	2.29	2.33	0.40	1.79	0.78	1.63
10	1.30	2.07	1.20	2.00	0.50	0.89
11	-2.18	2.73	1.33	1.73	1.90	1.97
12	-1.71	2.26	-1.00	1.96	0.07	1.07
13	0.25	2.54	-1.50	1.87	-0.25	1.53
14	-1.56	1.94	1.20	1.67	0.83	1.35

Pooled	70% - 80% SaO ₂	80% - 90% SaO ₂	90% - 100% SaO ₂
Mean Bias	0.16	0.21	0.21
Arms	2.00	1.87	1.29

Figure 1 Plot of difference (SpO₂ - SaO₂) versus artery blood gas (SaO₂) with linear regression fit and upper 95% and lower 95% limits of agreement of all subjects. Each color or symbol represents a different patient in the clinical study.



Difference plot of MbH TD-8255 Fingertip Pulse Oximeter and artery blood gas

GUIDANCE AND MANUFACTURER'S DECLARATION

Manufacturer's declaration-electromagnetic emissions		
The device is intended for use in the electromagnetic environment (for home and professional healthcare) specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment-guidance (for home and professional healthcare environment)
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

Manufacturer's declaration-electromagnetic immunity			
The device is intended for use in the electromagnetic environment (for home and professional healthcare) specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance (for home and professional healthcare environment)
Electrostatic discharge (ESD) IEC 61000-4-2	Contact: ±8 kV Air±2 kV, ±4 kV, ±8 kV, ±15 kV	Contact: ±8 kV Air±2 kV, ±4 kV, ±8 kV, ±15 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	Not applicable	Mains power quality should be that of a typical home and professional healthcare environment.
Surge IEC 61000-4-5	±0.5kV, ±1kV line(s) to line(s) ±0.5kV, ±1kV, ±2kV line(s) to earth	Not applicable	Mains power quality should be that of a typical home and professional healthcare environment.
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Voltage dips: 0% UT; 0.5 cycle 0% UT; 1 cycle 70% UT; 25/30 cycles Voltage interruptions: 0% UT; 250/300 cycle	Voltage dips: Not applicable Not applicable Not applicable Voltage interruptions: Not applicable	Mains power quality should be that of a typical home and professional healthcare environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50, 60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz and 60 Hz	The device power frequency magnetic fields should be at levels characteristic of a typical location in a typical home and professional healthcare environment.

NOTE UT is the a.c. mains voltage prior to application of the test level.

Manufacturer's declaration-electromagnetic immunity			
The device is intended for use in the electromagnetic environment (for home and professional healthcare) specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance (for home and professional healthcare environment)
Conducted RF IEC 61000-4-6	3 Vrms: 0.15 MHz – 80 MHz 6 Vrms: in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz	Not applicable	Portable and mobile RF communications equipment should be used no closer to any part of the device including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	10 V/m: 80 MHz – 2.7 GHz 80% AM at 1 kHz	10 V/m: 80 MHz – 2.7 GHz 80% AM at 1 kHz	Recommended separation distance: d = 1.2 √p d = 1.2 √p 80MHz to 800 MHz d = 2.3 √p 800MHz to 2,7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Recommended separation distance between portable and mobile RF communications equipment and the device			
The device is intended for use in an electromagnetic environment (for home and professional healthcare) in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz d = 1,2 √p	80 MHz to 800 MHz d = 1,2 √p	800 MHz to 2,7 GHz d = 2,3 √p
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Manufacturer's declaration-electromagnetic immunity							
Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment							
The device is intended for use in the electromagnetic environment (for home and professional healthcare) specified below. The customer or the user of the device should assure that it is used in such an environment.							
Test frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)	Compliance LEVEL (V/m) (for home and professional healthcare)
385	380 – 390	TETRA 400	Pulse modulation ^{b)} 18 Hz	1,8	0,3	27	27
450	430 – 470	GMRS 460, FRS 460	FM ^{c)} ±5 kHz deviation 1 kHz sine	2	0,3	28	28
710	704 – 787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9	9
745							
780							
810	800 – 960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	2	0,3	28	28
870							
930							
1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	2	0,3	28	28
1845							
1970							
2450							
2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0,3	28	28
5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9	9
5500							
5785							

NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.
^{a)} For some services, only the uplink frequencies are included.
^{b)} The carrier shall be modulated using a 50% duty cycle square wave signal.
^{c)} As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

Manufacturer's declaration-electromagnetic immunity				
Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields				
The device is intended for use in the electromagnetic environment (for home and professional healthcare) specified below. The customer or the user of the device should assure that it is used in such an environment.				
Frequencies	Test level [A/m]	Modulation	Dwell time [s]	Compliance LEVEL [A/m] (for home and professional healthcare)
30 kHz (a)	8	CW	3	8
134,2 kHz	65	Pulse modulation (b) 2,1 kHz	3	65 (c)
13,56 MHz	7,5	Pulse modulation (b) 2,1 kHz	3	7,5 (c)

Note:
(a) This test is applicable only to ME EQUIPMENT and ME SYSTEMS intended for use in the HOME AND PROFESSIONAL HEALTHCARE ENVIRONMENT.
(b) The carrier shall be modulated using a 50% duty cycle square wave signal.
(c) r.m.s, before modulation is applied.